



APPLICATION SHEET - Beckman Coulter® AU480, AU640, AU680, AU2700, AU5400, AU5800

MyCare™ Psychiatry Clozapine Assay

Reagent ID: CLZ-RGT

Clozapine

Configure assay parameters - Specific Test Parameters - General

General LIH ISE Range

Test Name: CLZ Type: Serum Operation: Yes

Sample:	Volume	<input type="text" value="4"/>	μL	Dilution	<input type="text" value="0"/>	μL	Pre-Dilution Rate:	<input type="text" value="1"/>
Reagents:	R1 Volume	<input type="text" value="100"/>	μL	Dilution	<input type="text" value="0"/>	μL	Min OD:	<input type="text" value=""/>
	R2 Volume	<input type="text" value="50"/>	μL	Dilution	<input type="text" value="0"/>	μL	Max OD:	<input type="text" value=""/>
Wavelength:	Pri.	<input type="text" value="600"/>		Sec.	<input type="text" value="(None)"/>		Reagent OD Limit:	
Method:		<input type="text" value="END1"/>					First L	<input type="text" value="-2"/>
Reaction slope:		<input type="text" value="+"/>					First H	<input type="text" value="3"/>
Measuring Point 1:	First	<input type="text" value="12"/>		Last	<input type="text" value="27"/>		Last L	<input type="text" value="-2"/>
Measuring Point 2:	First	<input type="text" value=""/>					Last H	<input type="text" value="3"/>
Linearity:		<input type="text" value=""/>	%				Dynamic Range:	
No Lag Time:		<input type="text" value=""/>					L	<input type="text" value="68"/>
							H	<input type="text" value="1500"/>
							Correlation Factor:	
							A	<input type="text" value="1"/>
							B	<input type="text" value="0"/>
							On-board stability period:	<input type="text" value="75"/>

Configure assay parameters - Specific Test Parameters - Range

General LIH ISE Range

Test Name: CLZ Type: Serum

Value/Flag:	<input type="text" value=""/>	Level L:	<input type="text" value=""/>	Level H:	<input type="text" value=""/>
Normal Ranges	Age L	Age H			
	Sex	Year	Month	Year	Month
1.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
2.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
3.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
4.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
5.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
6.	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>
7.	None Selected				
8.	Out of Range				
			L	H	
Panic Value:	<input type="text" value=""/>	<input type="text" value=""/>			
			Unit:	ng/mL	Decimal places:
					<input type="text" value="0"/>



APPLICATION SHEET - Beckman Coulter® AU480, AU640, AU680, AU2700, AU5400, AU5800

MyCare™ Psychiatry Clozapine Assay

Reagent ID: CLZ-RGT

Clozapine

Configure assay parameters - Calibration Parameters - Calibration Specific

General		ISE			
Test Name:	CLZ	Type:	Serum		
Calibration Type:	6AB	Formula:	EIA TYPE 1		
Counts:	2	Process:	Conc.		
Point 1:	Cal. No.	OD	CONC.*	Factor/OD-L	Factor/OD-H
Point 2:			0	-2.0000	3.0000
Point 3:			75	-2.0000	3.0000
Point 4:			121	-2.0000	3.0000
Point 5:			248	-2.0000	3.0000
Point 6:			539	-2.0000	3.0000
Point 7:			1423	-2.0000	3.0000
1 Point Cal. Point:		<input type="checkbox"/> with CONC-0	Advanced Calibration:		NO
MB Type Factor:		Calibration Stability Period:		30	

Configure assay parameters - QC Parameters - QC Specific - Check

Test Name	CLZ
Single Check Level	2SD
Multi Check Level	
QC Mode	Preset

Configure assay parameters - QC Parameters - QC Specific - Preset

Test Name	CLZ	Type:	Serum		
Level	Control**	Multi/Single	Mean	SD	Range
1	Low	Single	152	13	52
2	Med	Single	461	60	240
3	High	Single	907	118	472

____ User Defined

*Saladax Biomedical MyCare Psychiatry Calibrator Kit 2 Cat. No. MCP2-CAL; US Cat. No. MCP2-CAL-US

**Saladax Biomedical MyCare Psychiatry Control Kit 2 Cat. No. MCP2-CON; US Cat No. MCP2-CON-US

For further information refer to the package insert or contact Saladax Biomedical Customer Support at:

Telephone: +1 610 419-6731
Email: techsupport@saladax.com