Paclitaxel (MyPaclitaxel™) Assay

This package insert must be read carefully prior to product use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.


Customer Service:

Telephone: +1-610-419-6731

Fax: +1-484-547-0590

Key to Symbols Used

	<i>in vitro</i> Diagnostic Device		Consult Instructions for Use
	Batch Code (Lot)		Temperature Limitation
	Catalog Number		Use By
	Caution: Consult accompanying documents		Authorized Representative in the European Community
	Reagent 1		Manufacturer
	Reagent 2		Calibrator A

EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands

Saladax Biomedical, Inc.
 116 Research Drive
 Bethlehem, PA 18015 USA
www.saladax.com

Intended Use

The Saladax Paclitaxel (MyPaclitaxel) Assay is an *in vitro* diagnostic medical device intended for the quantitative determination of paclitaxel in human plasma using automated clinical chemistry analyzers as an aid in the management of paclitaxel therapy.

Summary and Explanation of the Test

Paclitaxel (Taxol®, Abraxane®) is a commonly used antineoplastic agent of the taxoid family with broad activity in patients with a variety of solid tumors such as ovarian cancer, breast cancer, non-small cell lung cancer (NSCLC), and AIDS-related Kaposi's Sarcoma^{1, 2}. Paclitaxel is not only used for advanced disease, but also in the adjuvant treatment of breast cancer³, and in early ovarian cancer⁴. As with most cytotoxic agents, paclitaxel has a narrow therapeutic window and inter-individual variability of pharmacokinetics (PK) is a major contributor to unpredictable toxicity.

Paclitaxel is typically administered at 135 - 175 mg/m² IV over 3 hours, every 3 weeks². However, this body surface area (BSA)-based dosing is associated with paclitaxel systemic exposure variability as high as 10-fold⁵⁻⁹. Genetic polymorphisms⁹⁻¹² and several other covariates (sex, BSA, and age)^{13, 14} have been shown to affect the elimination capacity of paclitaxel. Inter-patient variability in exposure could be reduced by using individual patient PK data to guide dosing. Several such studies using dose adjustment have been performed. The first one enrolled children with refractory acute leukemia¹⁵. Twelve hours into a 24-hour infusion, one arm had doses adjusted to achieve the target AUC of 2.9 mg•h/L. With adjustment, the resultant AUCs were between 75% and 125% of the target, whereas none of the unadjusted arm was in the target range ($p=0.021$). In another study, dose adjustments of paclitaxel were performed after first cycle administration of paclitaxel (175 mg/m² over 3 hours) for patients with advanced NSCLC to maintain a paclitaxel plasma concentration of 0.10 µmol/L for more than 15 hours¹⁶. Patients requiring adjustment to reach the target exposure duration decreased from 36% to 11% over 6 courses. The recently completed Central European Society for Anticancer Research (CESAR) Study of Paclitaxel Therapeutic Drug Monitoring (CEPAC-TDM) randomized 304 patients with advanced NSCLC receiving paclitaxel/carboplatin to either a standard dose (Arm A) or at a PK-guided dose (Arm B) for which subsequent PTX doses were adjusted to a $T_{C>0.05}$ of 26 - 31 hours based on the $T_{C>0.05}$ in the previous cycle. The proportion of patients with supratherapeutic PTX exposure was reduced from

41% in cycle 1 to 2% in cycle 6¹⁷. These studies demonstrate that paclitaxel therapeutic dose management (TDM) with individual dose adjustments is feasible, and may provide improved outcomes for patients undergoing chemotherapy with paclitaxel.

There is increasing evidence for the relationship between paclitaxel plasma concentrations and drug-associated toxicity and clinical outcome^{14, 17-20}. The major toxicity observed during paclitaxel treatment is bone marrow suppression (neutropenia, anemia, thrombocytopenia), followed by cumulative dose phenomenon such as neurotoxicity, with frequent peripheral neuropathy.

The PK parameter which has emerged as the most relevant predictor of hematological toxicity, neuropathy and also clinical outcome in some studies, is the time above a paclitaxel plasma concentration of 0.05 $\mu\text{mol/L}$ (less frequently 0.10 $\mu\text{mol/L}$ ¹⁸), or $T_{>0.05}$. The most frequent adverse effect was neutropenia^{7, 9, 12, 14, 19, 21, 22}. A large data-set from a multi-center pharmacology trial in 135 patients with ovarian cancer receiving paclitaxel/carboplatin treatment was published by Joerger *et al.* in 2007¹⁴. Paclitaxel $T_{>0.05}$ was significantly higher in patients with severe as compared to no or mild neutropenia. A paclitaxel $T_{>0.05}$ of between 24 and 30 hours was the best discriminator between mild and severe neutropenia. With respect to neuropathy, Mielke *et al.* showed a significant risk of peripheral neuropathy with paclitaxel $T_{>0.05} > 10.6$ hours ($p=0.023$)²³. In the larger randomized CEPAC-TDM study mentioned above PK-guided dosing resulted in a clinically relevant reduction of neuropathy compared to standard dosing (grade ≥ 2 neuropathy 14% vs. 27%, $p < 0.001$). Lower average PTX doses (132 vs. 197 mg/m^2 , $p < 0.001$) had no negative effect on clinical outcome in the overall study population - Arms A and B: objective response rate of 32% vs. 29% ($p=0.70$) and progression-free survival of 5.2 vs. 4.7 months (hazard ratio 1.1, 95% CI 0.8-1.4, $p=0.54$)¹⁷.

Paclitaxel time above threshold concentration was also related to clinical outcome, although to a lesser extent. Huizing *et al.* evaluated patients with advanced NSCLC, and noted different threshold values correlated with different outcomes. Those patients with longer paclitaxel $T_{>0.10}$ had a longer median survival as compared to patients with a shorter paclitaxel $T_{>0.10}$ (8.2 vs. 4.8 months, $p=0.06$), and paclitaxel $T_{>0.05}$ was correlated to a significantly improved time to disease progression (TTP) (82.0 vs 62.2 weeks, $p=0.05$)¹⁸. Similarly, Joerger *et al.* showed that in ovarian cancer patients receiving 3-weekly paclitaxel/carboplatin treatments, $T_{>0.05}$ correlated to improved TTP (85.2 versus 63.3 weeks; $p=0.05$)¹⁴.

The goals of paclitaxel TDM are to identify those patients that are at risk for toxicity. This information can be used 1) to make a decision whether to prophylactically administer granulocyte colony-stimulating factor (GCSF) to prevent neutropenia and/or 2) to *guide* dose administration to an effective exposure, while remaining below plasma concentrations that give rise to toxicity.

Currently PK-guided TDM for paclitaxel is not a clinical reality in part due to the lack of rapid and cost-effective methods to monitor cancer drug plasma concentrations. Monitoring *is possible* by chromatographic methods such as liquid chromatography mass spectroscopy (LC-MS), but the cost per test for these methodologies is high and requires extensive sample pretreatment. The sample pretreatment steps (extraction into organic solvent, drying, and resuspension into aqueous solvent), can introduce significant error in the final measurement. By comparison, the MyPaclitaxel Assay provides a convenient, cost-effective, and timely tool to aid oncologists in paclitaxel dose management.

Principles of the Procedure

The MyPaclitaxel Assay (US Patent No. 7,175,993) is a homogeneous two-reagent nanoparticle agglutination immunoassay used to measure the concentration of paclitaxel in human plasma. It is based upon the principle of measuring changes in scattered light, or apparent absorbance, which result when antibody coated nanoparticles aggregate. This agglutination is measured at wavelengths between 400 and 650 nm using automated clinical analyzers. In this technology, multivalent drug-conjugates serve as binding partners to antibodies selective for paclitaxel, which are covalently attached to the surface of nanoparticles. In the absence of free paclitaxel, this reaction creates large aggregates, resulting in a solution that scatters incident light and leads to an increase in the observed absorption of the solution. When a sample containing paclitaxel is introduced, the agglutination reaction is partially inhibited. Antibody bound to drug in the sample is no longer available to promote nanoparticle aggregation, resulting in less scattering of incident light and lower observed absorption of the solution. This results in an inhibition curve with increasing paclitaxel: no drug in the sample results in maximum absorbance and high drug levels in the sample result in a minimal amount of absorbance. Monitoring the changes in scattered light or absorbance as a function of drug levels results in a concentration-dependent curve^{24, 25}.

Reagents

MyPaclitaxel Assay (REF PTX-RGT) includes:

1. **Reaction Buffer** (R1): contains agglutination enhancers, protein, and buffer.
2. **Nanoparticle Reagent** (R2): contains monoclonal antibody bound to nanoparticles in a buffered solution.

Additional Materials Required But Not Provided:

(REF) PTX-CAL – MyPaclitaxel Calibrator Kit

(REF) PTX-CON – MyPaclitaxel Control Kit

(REF) PTX-CON-C111 – MyPaclitaxel Control Kit for cobas® c111

Materials Provided But Not Required:

(REF) GEN-DEC – MyCare Drug Exposure Calculator

Precautions and Warnings

For In Vitro Diagnostic Use Only.

Exercise the normal precautions required for handling all laboratory reagents.

All components of the MyPaclitaxel Assay contain less than 0.1% sodium azide. For the specific listing, refer to the reagent section of this package insert. Avoid contact with skin and mucous membranes. Flush affected areas with copious amounts of water. Seek immediate medical attention if reagents are ingested or come into contact with eyes. When disposing of such reagents, always flush with large amounts of water to prevent accumulation of azide.

Handling and Storage Instructions

Store reagents, calibrators and controls refrigerated at 2-8°C (35-46°F). Do not freeze.

Before use, mix the Nanoparticle Reagent (R2) by gently inverting the R2 reagent vial three to five times, avoiding the formation of bubbles.

Indications of Stability

Reagents are stable until the expiration date when stored and handled as directed.

Sample Collection and Handling

Plasma (EDTA) specimens may be used with the MyPaclitaxel Assay. Collect a minimum of 2 mL of blood into an EDTA tube. Collect the blood sample by venipuncture or through a peripheral IV line. Draw the sample **24 ± 6 (18-30) hours after the start of infusion**. The start time of the infusion and actual sampling times should be recorded.

Collected samples should be placed on ice immediately and be processed to plasma within 12 hours of collection. Draw off plasma from the top of the tube, avoiding the cell layer (contamination of plasma with blood cells may interfere with results), and transfer to a capped secondary tube. Store the samples at 2-8°C for immediate testing (within 7 days) or freeze ($\leq -20^{\circ}\text{C}$) if storing for a longer period. If the samples will be shipped to an outside testing laboratory, they may be stored at 2-8°C or frozen and shipped on ice packs for delivery within 7 days.

Assay Procedure

Refer to the instrument specific Application Sheet for instructions and parameters before performing the assay.

Specimen Dilution Procedure

Samples containing paclitaxel concentrations greater than 320 ng/mL can be diluted 1:2 for a reportable assay range up to 960 ng/mL. Specimens out of range can be manually diluted 1:2 with paclitaxel Calibrator A (CALA ; one part specimen and two parts Calibrator A) prior to re-testing.

Calibration

The MyPaclitaxel Assay produces a calibration curve with a 0 to 320 ng/mL range using the MyPaclitaxel Calibrator Kit. The minimum detectable concentration of paclitaxel in plasma for the MyPaclitaxel Assay is 17 ng/mL.

Verify the assay calibration by testing MyPaclitaxel Controls.

Calibration Frequency

Calibration is recommended:

- After a lot change of calibrator or reagent (kit),
- After performing major instrument maintenance or repairs,
- As required by each laboratory's quality control procedures.

Quality Control

The MyPaclitaxel Control Kit contains three levels of controls with low, medium and high paclitaxel concentrations.

Each laboratory should establish its own control ranges and frequency of quality control. Good laboratory practice suggests that at least two control levels be tested each day patient samples are assayed and each time a calibration is performed. Reassess control targets and ranges following a lot change of reagent (kit), calibrator, or control.

Results and Expected Values

The analyzer software calculates a best fit non-linear curve equation that is used to generate a calibration curve that ranges from 0 to 320 ng/mL paclitaxel. Concentrations of drug in unknown samples are calculated from this calibration curve using absorbance values generated for each sample.

Limitations of the Procedure

As with all analyte determinations, the MyPaclitaxel Assay should be used in conjunction with information available from clinical evaluation and other diagnostic procedures.

Performance characteristics for the MyPaclitaxel Assay have not been established for body fluids other than human plasma containing EDTA.

No significant interferences were observed from samples with the following conditions:

	Level	
Rheumatoid Factor	500 IU/mL	
Total Protein Matrix Effect	15.0 g/dL	150 g/L
Icteric Interference	96 mg/dL	1,642 μ mol/L
Lipemic Interference	603 mg/dL	6.6 mmol/L
Hemolysate	1,000 mg/dL	

Paclitaxel is not stable in whole blood. Hemolysis should be avoided.

Paclitaxel metabolites 6 α -hydroxy-paclitaxel and 3'*p*-hydroxy-paclitaxel have 3.4% and 3.0% cross-reactivity, respectively, in the MyPaclitaxel Assay. These metabolites are rapidly broken down further in the liver, and will not accrue to detectable quantities in the blood at the recommended sampling timepoint. Maximum metabolite concentrations in human blood are < 10% of the total paclitaxel for 6 α -hydroxy-paclitaxel and < 5% for 3'*p*-hydroxy-paclitaxel, which would contribute < 0.34% and < 0.15% positive bias, respectively.

Paclitaxel will epimerize under alkaline conditions to 7-*epi*-paclitaxel. Patient plasma samples do not contain the epimer upon collection, but prolonged storage at 2-8°C (> 7 days) or room temperature (> 6 hours) will result in epimer formation. Proper handling of patient samples is essential to avoid biased results.

As with any assay using mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) in the sample. Samples containing such antibodies can potentially produce erroneous results, which are inconsistent with the patient's clinical profile. If HAMA interference is suspected, please contact Saladax Technical Service for assistance.

Expected Values

No precise relationship between paclitaxel plasma concentrations and antineoplastic efficacy has been established. A number of reports have demonstrated wide pharmacokinetic variability of the drug and some reports in the literature suggest that a time above a threshold plasma concentration of 0.05 μ mol/L ($T_{C>0.05}$: 42.7 ng/mL) increases the risk of toxicity. In the literature reported $T_{C>0.05}$ values range from 5-55 hours^{20, 22, 23}, and plasma concentrations above 0.05 μ mol/L for periods greater than 30 hours have been associated with

increased risk of toxicity in the form of grade 3/4 neutropenia^{19, 20}. The $T_{C>0.05}$ can be calculated using commercial software for the determination of pharmacokinetic parameters or with the MyCare Drug Exposure Calculator (GEN-DEC). Blood samples to determine $T_{C>0.05}$ can be drawn using existing limited sampling strategies which involve 2 draw times^{26, 27}. However, Joerger *et al.* have developed a limited sampling strategy using only one draw^{28, 29}. From a data set of 135 patients with draw times of 18-30 hours¹⁴, the measured plasma concentrations were below 320 ng/mL, with the majority between 40 and 80 ng/mL, and only 12% below 20 ng/mL.

Paclitaxel drug concentrations should not be the only means of therapeutic drug management. The assay should be used in conjunction with information available from clinical evaluations and other diagnostic procedures. The patient's current and past medical condition, the complexity of the clinical state, individual differences in sensitivity to paclitaxel and toxic effects of paclitaxel, co-administration of other drugs, and a number of other factors may result in different optimal blood concentrations of paclitaxel for any individual. Each patient should undergo comprehensive clinical assessment prior to modification of the treatment plan and clinicians should carefully monitor patients during therapy initiation and dose adjustments. Given the heterogeneity of the patient's clinical state, clinicians should establish a desired therapeutic management range based on their own experience as well as each patient's clinical requirements. Paclitaxel concentrations for individual patients should be determined using a single, consistent method to minimize confounding effects associated with cross-reactivity and recognition of metabolites.

Specific Performance Characteristics

Typical performance data for the MyPaclitaxel Assay obtained on a Beckman (formerly Olympus) AU400™ or Beckman Coulter AU480™ are shown below. Results obtained in individual laboratories may differ from these data.

Precision

Precision was determined as described in CLSI Guideline EP5-A2.

Low, medium and high paclitaxel controls and 4 patient sample pools containing varying paclitaxel concentrations were assayed in duplicate twice a day for 20 days at one site using 2 assay lots and two analyzers. The means were

determined and SD and % CVs were calculated for repeatability and within-laboratory precision.

The following are representative results from two assay lots at one site.

Sample Type	Assigned Value (ng/mL)	N	Mean (ng/mL)	Repeatability		Within-laboratory precision	
				SD	CV	SD	CV
Controls	60	320	61	2.0	3.3%	3.9	6.4%
	120	336	124	2.4	2.0%	4.8	3.9%
	200	336	202	3.1	1.5%	6.0	3.0%
Human Plasma	50	336	50	1.7	3.4%	5.0	9.9%
	100	336	98	2.0	2.0%	6.0	6.1%
	180	336	177	3.6	2.0%	8.3	4.7%
	280	336	278	4.1	1.5%	9.6	3.5%

Limit of Quantitation (LoQ)

LoQ was determined as described in CLSI Guideline EP17-A.

The LoQ of an assay is the lowest analyte concentration that can be measured with acceptable accuracy and precision. For the MyPaclitaxel Assay, the LoQ was defined as the lowest paclitaxel concentration at which the total analytical error was $\leq 35\%$. The LoQ determined was 24 ng/mL.

Limit of Detection (LoD)

LoD was determined as described in CLSI Guideline EP17-A.

The LoD of an assay is the lowest analyte concentration that can be detected with stated confidence, but not necessarily measured quantitatively. For the MyPaclitaxel assay, the LoD was defined as the lowest paclitaxel concentration at which 95% of the results exceeded the assay Limit of Blank (LoB), where the assay LoB was the 95th percentile of the results generated from the analysis of specimens that did not contain paclitaxel. The LoD determined was 17 ng/mL.

Specificity

Paclitaxel metabolites and Structurally Related Compounds

Indicated concentrations of each of the following paclitaxel metabolites or structurally related compounds were added to 150 ng/mL paclitaxel-spiked plasma and were assayed using the MyPaclitaxel Assay. The cross-reactivity in the assay for each compound is given below.

Compound	Cross-Reactivity	
7-epi-paclitaxel	100 ng/mL	200 %
6 α -hydroxy-paclitaxel	500 ng/mL	3.4 %
3' <i>p</i> -hydroxy-paclitaxel	500 ng/mL	3.0 %
Baccatin III	10,000 ng/mL	< 0.1 %
10-deacetyl-paclitaxel	100 ng/mL	113 %
10-deacetyl-baccatin	10,000 ng/mL	1.2 %

Common Co-Administered Drugs

Unless otherwise indicated, 100,000 ng/mL of each compound was spiked into paclitaxel free plasma or plasma spiked with 150 ng/mL Paclitaxel. The cross-reactivity of all compounds tested in the assay was $\leq 0.1\%$.

1-methylxanthine	Cromolyn	Phenobarbital**
3-methylxanthine	Cyclophosphamide	Phenytoin
5-FU	Dexamethasone	Prednisolone
5-methyl-1H-benzimidazole	Doxorubicin	Prednisone
6-mercaptopurine	Erythromycin	Procainamide
7-methylxanthine*	Folic acid	Prochlorperazine
8-chlorotheophylline	Formoterol	Quinidine sulfate
Acetaminophen	Gentamycin	Retinol
Acetylsalicylic acid	Heparin Lithium Salt	Rifampicin
Adenine	Heparin Sodium Salt	Salicylic acid
Albuterol (Salbutamol)	Hydrocortisone	Spectinomycin
Allopurinol	Hypoxanthine	Streptomycin
Amikacin	Ibuprofen	Theobromine
Amoxicillin	Ipratropium	Theophylline
Ampicillin	Irinotecan	Tobramycin
Azathioprine	Kanamycin A	Valproic acid
Benzimidazole	Leucovorin	Vancomycin

Caffeine	Lidocaine	Vitamin B1
Calcium carbonate	Methotrexate	Vitamin B2
Calcium citrate	Morphine**	Vitamin B3
Capecitabine	N-acetylprocainamide	Vitamin B6
Carbamazepine	Nicotine	Vitamin B12
Carbamazepine-10,11-epoxide*	Nornicotine	Vitamin C
Carboplatin	Oxaliplatin	Vitamin D3
Ceftriaxone	Oxypurinol	Vitamin E
Cephalosporin	Penicillin G	Vitamin K1
Cisplatin	Pentobarbital**	Xanthine
Cotinine	Pentoxifylline	β-carotene

*50,000 ng/mL **10,000 ng/mL

Recovery

Assay recovery is the degree to which an analyte result reported by an assay agrees with the true analyte concentration in the specimen. For the MyPaclitaxel Assay recovery was defined as the percent recovery of paclitaxel in normal plasma specimens spiked with paclitaxel relative to the theoretical paclitaxel concentration.

To assess recovery, paclitaxel was spiked at 4 concentrations: 50, 70, 140, and 280 ng/mL into 4 normal K₂EDTA plasmas. The percent deviations at each spike level among the calibrators lots and analyzers evaluated ranged from 87%-113%.


Linearity

Assay linearity was determined as described in CLSI Guideline EP6-A, which states that quantitative results in the linear range of the assay are directly proportional to the concentration of the analyte in the sample. In the linear range of the MyPaclitaxel at levels >30 ng/mL the deviation from linearity was within ± 2%, individual replicate deviation was <14% and mean deviation was less than 7%. At 30 ng/mL the the greatest deviation from linearity was 15%, the greatest individual deviation for linearity samples was less than 23%, and the greatest mean deviation for linearity samples was less than 17%.

To determine assay linearity, 11 samples were prepared according to CLSI EP6-A. The assay was linear from 20-330 ng/mL.

Method Comparison

A comparison between the Saladax MyPaclitaxel Assay and an HPLC-UV method was performed using 57 human plasma samples obtained from patients receiving paclitaxel therapy. The range of paclitaxel concentrations in the MyPaclitaxel Assay was 28 – 257 ng/mL with a mean of 112 ng/mL. The range of paclitaxel concentrations for the HPLC-UV reference method was 20 – 234 ng/mL with a mean of 111 ng/mL. Results of the Deming regression analysis are below.

Package Insert References:

1. Abraxane, *Abraxane for Injectable Suspension (paclitaxel protein-bound particles for injectable suspension) (albumin-bound)*, in *Package Insert*, 2007, Abraxis: Los Angeles (CA).
2. Taxol, *Taxol (Paclitaxel) Injection*, in *Package Insert*, 2011, Bristol-Myers Squibb: Princeton (NJ).
3. Sparano, J.A., et al., *Weekly paclitaxel in the adjuvant treatment of breast cancer*. New England Journal of Medicine, 2008. **358**(16): p. 1663-71.
4. Ozols, R.F., et al., *Phase III trial of carboplatin and paclitaxel compared with cisplatin and paclitaxel in patients with optimally resected stage III ovarian cancer: a Gynecologic Oncology Group study*. Journal of Clinical Oncology, 2003. **21**(17): p. 3194-3200.
5. Sonnichsen, D.S., et al., *Saturable pharmacokinetics and paclitaxel pharmacodynamics in children with solid tumors*. J Clin Oncol, 1994. **12**(3): p. 532-8.
6. Henningsson, A., et al., *Mechanism-based pharmacokinetic model for paclitaxel*. J Clin Oncol, 2001. **19**(20): p. 4065-73.
7. Mould, D.R., et al., *Population analysis of a 24-h paclitaxel infusion in advanced endometrial cancer: a gynaecological oncology group study*. Br J Clin Pharmacol, 2006. **62**(1): p. 56-70.
8. Huizing, M.T., et al., *Pharmacokinetics of paclitaxel and metabolites in a randomized comparative study in platinum-pretreated ovarian cancer patients*. J Clin Oncol, 1993. **11**(11): p. 2127-35.
9. Jiko, M., et al., *Pharmacokinetics and pharmacodynamics of paclitaxel with carboplatin or gemcitabine, and effects of CYP3A5 and MDR1 polymorphisms in patients with urogenital cancers*. Int J Clin Oncol, 2007. **12**(4): p. 284-90.
10. Sissung, T.M., et al., *Association of ABCB1 genotypes with paclitaxel-mediated peripheral neuropathy and neutropenia*. Eur J Cancer, 2006. **42**(17): p. 2893-6.
11. Henningsson, A., et al., *Association of CYP2C8, CYP3A4, CYP3A5, and ABCB1 polymorphisms with the pharmacokinetics of paclitaxel*. Clin Cancer Res, 2005. **11**(22): p. 8097-104.
12. Nakajima, M., et al., *Pharmacokinetics of paclitaxel in ovarian cancer patients and genetic polymorphisms of CYP2C8, CYP3A4, and MDR1*. J Clin Pharmacol, 2005. **45**(6): p. 674-82.
13. Joerger, M., et al., *Quantitative effect of gender, age, liver function, and body size on the population pharmacokinetics of Paclitaxel in patients with solid tumors*. Clin Cancer Res, 2006. **12**(7 Pt 1): p. 2150-7.
14. Joerger, M., et al., *Population pharmacokinetics and pharmacodynamics of paclitaxel and carboplatin in ovarian cancer patients: a study by the European organization for research and treatment of cancer-pharmacology and molecular mechanisms group and new drug development group*. Clin Cancer Res, 2007. **13**(21): p. 6410-8.
15. Woo, M.H., et al., *Phase I targeted systemic exposure study of paclitaxel in children with refractory acute leukemias*. Clinical Cancer Research, 1999. **5**(3): p. 543-549.

16. de Jonge, M.E., et al., *Bayesian pharmacokinetically guided dosing of paclitaxel in patients with non-small cell lung cancer*. Clinical Cancer Research, 2004. **10**(7): p. 2237-2244.
17. Joerger, M., et al., *Open-label, randomized study of individualized, pharmacokinetically (PK)-guided dosing of paclitaxel combined with carboplatin in advanced non-small cell lung cancer (NSCLC) patient*. J Clin Oncol, 2015. **33**(suppl; abstr 8051).
18. Huizing, M.T., et al., *Pharmacokinetics of paclitaxel and carboplatin in a dose-escalating and dose-sequencing study in patients with non-small-cell lung cancer*. Journal of Clinical Oncology, 1997. **15**(1): p. 317-329.
19. Miller, A.A., et al., *Prospective evaluation of body surface area as a determinant of paclitaxel pharmacokinetics and pharmacodynamics in women with solid tumors: Cancer and Leukemia Group B Study 9763*. Clinical Cancer Research, 2004. **10**(24): p. 8325-8331.
20. Ohtsu, T., et al., *Clinical pharmacokinetics and pharmacodynamics of paclitaxel: A 3-hour infusion versus a 24-hour infusion*. Clinical Cancer Research, 1995. **1**(6): p. 599-606.
21. Kobayashi, M., et al., *Pharmacokinetic study of weekly administration dose of paclitaxel in patients with advanced or recurrent gastric cancer in Japan*. Gastric Cancer, 2007. **10**(1): p. 52-7.
22. Gianni, L., et al., *Nonlinear pharmacokinetics and metabolism of paclitaxel and its pharmacokinetic/pharmacodynamic relationships in humans*. Journal of Clinical Oncology, 1995. **13**(1): p. 180-190.
23. Mielke, S., et al., *Association of Paclitaxel pharmacokinetics with the development of peripheral neuropathy in patients with advanced cancer*. Clin Cancer Res, 2005. **11**(13): p. 4843-50.
24. Li, Z., et al., *New synthesis and characterization of (+)-lysergic acid diethylamide (LSD) derivatives and the development of a microparticle-based immunoassay for the detection of LSD and its metabolites*. Bioconjug Chem, 1997. **8**(6): p. 896-905.
25. McNally, A.J., et al., *An online immunoassay for LSD: comparison with GC-MS and the Abuscreen RIA*. J Anal Toxicol, 1996. **20**(6): p. 404-8.
26. Huizing, M.T., et al., *Limited sampling strategies for investigating paclitaxel pharmacokinetics in patients receiving 175 mg/m² as a 3-hour infusion*. Clinical Drug Investigation, 1995. **9**(6): p. 344-353.
27. Stroud, J.R., P. Muller, and G.L. Rosner, *Optimal sampling times in population pharmacokinetic studies*. Journal of the Royal Statistical Society Series C-Applied Statistics, 2001. **50**(part 3): p. 345-359.
28. Joerger, M., et al., *Evaluation of a pharmacology-driven dosing algorithm of 3-weekly paclitaxel using therapeutic drug monitoring: a pharmacokinetic-pharmacodynamic simulation study*. Clin Pharmacokinet, 2012. **51**(9): p. 607-17.
29. CESAR-EWIV, *An Open-Label, Randomized, Parallel Group Study of Patients Treated with Paclitaxel with Standard Dosing versus Pharmacokinetic Guided Dose Adjustment in Patients with Advanced NSCLC*. in CESAR Study of Paclitaxel Therapeutic Drug Monitoring: CEPAC-TDM, clinical study. CESAR C-III-002 / EUDRACT 2010-023688-16, Completed.

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