## MyCare Insite Clozapine Test

### INDIKATIONER FOR BRUG

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

### REFERENCER

1. Novartis Pharmaceuticals Corporation. Clozaril (clozapine) prescribing information. 2015
2. Hasegawa M, Gutierrez-Esteinou R, Way L, Meltzer HY. Relationship between clinical efficacy and clozapine concentrations in plasma in schizophrenia: effect of smoking. J Clin Psychopharmacol. 1993;13(6):383-390.
3. Kronig MH, Munne RA, Szymanski S, et al. Plasma clozapine levels and clinical response for treatment-refractory schizophrenic patients. Am J Psychiatry. 1995;152(2):179-182.
4. Miller DD, Fleming F, Holman TL, Perry PJ. Plasma clozapine concentrations as a predictor of clinical response: a follow-up study. J Clin Psychiatry. 1994;55 Suppl B:117-121.

5. Perry PJ, Miller DD, Arndt SV, Cadoret RJ. Clozapine and norclozapine plasma concentrations and clinical response of treatment-refractory schizophrenic patients. *Am J Psychiatry*. 1991;148(2):231-235.
6. Potkin SG, Bera R, Gulasekaram B, et al. Plasma clozapine concentrations predict clinical response in treatment-resistant schizophrenia. *J Clin Psychiatry*. 1994;55:133-136.
7. Spina E, Avenoso A, Facciola G, et al. Relationship between plasma concentrations of clozapine and norclozapine and therapeutic response in patients with schizophrenia resistant to conventional neuroleptics. *Psychopharmacology (Berl)*. 2000;148(1):83-89.
8. VanderZwaag C, McGee M, McEvoy JP, Freudenreich O, Wilson WH, Cooper TB. Response of patients with treatment-refractory schizophrenia to clozapine within three serum level ranges. *American Journal of Psychiatry*. 1996;153(12):1579-1584.
9. Lehman AF, Lieberman JA, Dixon LB, et al. Practice guideline for the treatment of patients with schizophrenia, second edition. *Am J Psychiatry*. 2004;161(2 Suppl):1-56.
10. Taylor D, Barnes TE, Young A. *The Maudsley Prescribing Guidelines in Psychiatry*. 13th ed. Hoboken, NJ: Wiley; 2018.
11. The American Psychiatric Association Practice Guideline For The Treatment Of Patients With Schizophrenia. 3rd ed. [ebook] Washington, DC: American Psychiatric Association Publishing., p.90. Available at: <https://psychiatryonline.org/doi/book/10.1176/appi.books.97808904248> [Accessed 8 September 2020]
12. Meyer JM, Stahl SM. *The Clozapine Handbook: Stahl's Handbooks*. Cambridge: Cambridge University Press; 2019.
13. Freudenreich O, McEvoy JP. Guidelines for prescribing clozapine in schizophrenia. In: Post TW, ed. *UpToDate*. Waltham, MA.: UpToDate; 2019. [www.uptodate.com](http://www.uptodate.com). Accessed February 21, 2019.
14. Buchanan RW, Kreyenbuhl J, Kelly DL, et al. The 2009 schizophrenia PORT psychopharmacological treatment recommendations and summary statements. *Schizophr Bull*. 2010;36(1):71-93.
15. Hiemke C, Bergemann N, Clement HW, et al. Consensus Guidelines for Therapeutic Drug Monitoring in Neuropsychopharmacology: Update 2017. *Pharmacopsychiatry*. 2018;51(1-02):9-62.
16. Schoretsanitis G, Kane JM, Correll CU, et al. Blood Levels to Optimize Antipsychotic Treatment in Clinical Practice: A Joint Consensus Statement of the American Society of Clinical Psychopharmacology and the Therapeutic Drug Monitoring Task Force of the Arbeitsgemeinschaft für Neuropsychopharmakologie und Pharmakopsychiatrie. *J Clin Psychiatry*. 2020;81(3)
17. Velligan DI, Weiden PJ, Sajatovic M, et al. Strategies for addressing adherence problems in patients with serious and persistent mental illness: recommendations from the expert consensus guidelines. *J Psychiatr Pract*. 2010;16(5):306-324.
18. Higashi K, Medic G, Littlewood KJ, Diez T, Granstrom O, De Hert M. Medication adherence in schizophrenia: factors influencing adherence and consequences of nonadherence, a systematic literature review. *Ther Adv Psychopharmacol*. 2013;3(4):200-218.
19. Velligan D, Weiden PJ, Sajatovic M, et al. Assessment of Adherence Problems in Patients with Serious and Persistent Mental Illness: Recommendations from the Expert Consensus Guidelines. *Journal of Psychiatric Practice*. 2010;16(1):34-45.
20. Mennicken S, Sobarzo A, Vega M, et al. Determination of clozapine in serum of patients with schizophrenia as a measurement of medication compliance. *Int J Psychiatry Clin Pract*. 2010;14(1):41-46.
21. Rostami-Hodjegan A, Amin AM, Spencer EP, Lennard MS, Tucker GT, Flanagan RJ. Influence of dose, cigarette smoking, age, sex, and metabolic activity on plasma clozapine concentrations: a predictive model and nomograms to aid clozapine dose adjustment and to assess compliance in individual patients. *J Clin Psychopharmacol*. 2004;24(1):70-78.
22. National Council Medical Director Institute (2018). *Medication Matters Causes and Solutions to Medication Non-Adherence*. Washington, DC: National Council for Behavioral Health.
23. Jakobsen MI, Larsen JR, Svensson CK, et al. The significance of sampling time in therapeutic drug monitoring of clozapine. *Acta Psychiatr Scand*. 2016;6(10):12673.
24. WHO Guidelines on Drawing Blood: Best Practices in Phlebotomy. Geneva: World Health Organization; 2010. 7, Capillary sampling. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK138654/>
25. Clinical Laboratory Standards Institute (CLSI). *Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard Sixth Edition*. CLSI document GP42-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
26. Gaertner I, Gaertner HJ, Vonthein R, Dietz K. Therapeutic drug monitoring of clozapine in relapse prevention: a five-year prospective study. *J Clin Psychopharmacol*. 2001;21(3):305-310. CLSI
27. CLSI. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition*. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute, 2014
28. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012
29. CLSI. *Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline – 2<sup>nd</sup> Edition*. CLSI document EP06. Wayne, PA: Clinical and Laboratory Standards Institute; 2020.
30. Saladax Biomedical, Inc. *MyCare Psychiatry Clozapine Assay Kit [package insert]* 2019.
31. CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples*. 3<sup>rd</sup> ed. CLSI guideline EP09c. Wayne, PA: Clinical and Laboratory Standards Institute, 2018.

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© 2024 Saladax Biomedical, Inc.