My5-FU

5-Fluorouracil (My5-FU™) Assay

Kundtjänst:

 Telefon: +1 (610) 419-6731
 Fax: +1 (484) 547-0590

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
Lutning = 1,017
y-intercept = 6,78
Korrelationskoefficient (R) = 0,9925

Referenser i produktbladet:

1. Teva Parenteral Medicines, Inc. Aduracil. Full Prescribing Information. 2015
2. Pfister DG, Spencer S, Brizel DM, et al. Head and neck cancers, version 1.2015. *J Natl Compr Canc Netw*. 2015;13:847-55.
3. Meyerhardt JA, Mayer RJ. Systemic therapy for colorectal cancer. *N Engl J Med*. 2005;352:476-487.
4. Gamelin E, Boisdron-Celle M. Dose monitoring of 5-fluorouracil in patients with colorectal or head and neck cancer – status of the art. *Crit Rev Oncol Hematol*. 1999;30:71-79.
5. Gamelin E, Boisdron-Celle M, Delva R, et al. Long-term weekly treatment of colorectal metastatic cancer with fluorouracil and leucovorin: results of a multicentric prospective trial of fluorouracil dosage optimization by pharmacokinetic monitoring in 152 patients. *J Clin Oncol*. 1998;16(4):1470-8.
6. Santini J, Milano G, Thyss A, et al. 5-FU therapeutic monitoring with dose adjustment leads to an improved therapeutic index in head and neck cancer. *Br J Cancer*. 1989;59(2):287-90.
7. Milano G, Etienne MC, Renee N, Thyss A, Schneider M, Ramaioli A, Demard F. et al. Relationship between fluorouracil systemic exposure and tumor response and patient survival. *J Clin Oncol*. 1994;12(6):1291-5.
8. Fety R, Etienne MC, Renee N, et al. Clinical impact of pharmacokinetically-guided dose adaptation of 5-fluorouracil: results from a multicentric randomized trial in patients with locally advanced head and neck carcinomas. *Clin Cancer Res*. 1998;4:2039-2045.
9. Levy E, Piedbois P, Buyse M, et al. Toxicity of fluorouracil in patients with advanced colorectal cancer: Effect of administration schedule and prognostic factors. *J Clin Oncol*. 1998;16:3537–3541.
10. Milano G, McLeod HL. Can dihydropyrimidine dehydrogenase impact 5-fluorouracil-based treatment? *Eur J Cancer*. 2000;36:37-42.
11. Jansman FG, Sleijfer DT, Coenen JL, De Graaf JC, Brouwers JR. Risk factors determining chemotherapeutic toxicity in patients with advanced colorectal cancer. *Drug Saf*. 2000;23:255-78.
12. Diasio RB. Sorivudine and 5-fluorouracil; a clinically significant drug-drug interaction due to inhibition of dihydropyrimidine dehydrogenase. *Br J Clin Pharmacol*. 1998;46:1-4.
13. Johnson MR, Diasio RB. Importance of dihydropyrimidine dehydrogenase (DPD) deficiency in patients exhibiting toxicity following treatment with 5-fluorouracil. *Adv Enzyme Regul*. 2001;41:151-7.
14. Ezzeldin H, Diasio R. Dihydropyrimidine Dehydrogenase Deficiency, a Pharmacogenetic Syndrome Associated with Potentially Life-Threatening Toxicity Following 5-Fluorouracil Administration. *Clin Colorectal Cancer*. 2004;4:181-189.
15. Milano G, Etienne M-C. Individualizing Therapy with 5-Fluorouracil Related to Dihydropyrimidine Dehydrogenase: Theory and Limits. *Ther Drug Monit*. 1996;18:335-340.
16. Capitain O, Asevoaia A, Boisdron-Celle M, Poirier AL, Morel A, Gamelin E. Individual fluorouracil dose adjustment in folfox based on pharmacokinetic follow-up compared with conventional body-area-surface dosing: A phase II, proof-of-concept study. *Clin Colorectal Cancer*. 2012;11:263-267.
17. Gamelin E, Jacob J, Merrouche Y, et al. Individual 5-fluorouracil dose adjustment based on pharmacokinetic follow-up compared with conventional dosage: results of a multicenter phase III randomized trial in patients with metastatic colorectal cancer. *J Clin Oncol*. 2008;26:2099-2105.
18. Kaldate RR, Haregewoin A, Grier CE, Hamilton SA, McLeod HL. Modeling the 5-fluorouracil area under the curve versus dose relationship to develop a pharmacokinetic dosing algorithm for colorectal cancer patients receiving FOLFOX6. *Oncologist*. 2012;17:296-302.

19. Kline CL, Schiccitano A, Zhu J, et al. Personalized dosing via pharmacokinetic monitoring of 5-fluorouracil might reduce toxicity in early- or late-stage colorectal cancer patients treated with infusional 5-fluorouracil-based chemotherapy regimens. *Clin Colorectal Cancer* 2014;13:119-126.
20. Patel JN, O'Neil BH, Deal AM, et al. A community-based multicenter trial of pharmacokinetically guided 5-fluorouracil dosing for personalized colorectal cancer therapy. *Oncologist*. 2014;19:959-965.
21. Braiteh FS, Salamone SJ, Li Y, et al. Pharmacokinetic (PK)-guided optimization of 5-fluorouracil (5FU) exposure in colorectal cancer (CRC) patients: U.S.-based clinical practices experience. *J Clin Oncol*. 2014;32:(suppl: abstr 3574).
22. Wilhelm M, Mueller L, Miller MC, et al. Prospective, multi-center study of 5-fluorouracil therapeutic drug management in metastatic colorectal cancer treated in routine clinical practice. *Clinical Colorectal Cancer*. 2016;15:381-388.
23. Saam J, Critchfield GC, Hamilton SA, Roa BB, Wenstrup RJ, Kaldate RR. Body surface area-based dosing of 5-fluorouracil results in extensive interindividual variability in 5-fluorouracil exposure in colorectal cancer patients on FOLFOX regimens. *Clin Colorectal Cancer*. 2011;10:203-206.
24. Hillcoat BL, McCulloch PB, Figueredo AT, et al. Clinical response and plasma levels of 5-fluorouracil in patients with colonic cancer treated by drug infusion. *Br J Cancer*. 1978;38:719-724.
25. Seitz JF, Cano JP, Rigault JP, et al. Chimiothérapie des cancers digestifs étendus par le 5-Fluorouracile: relations entre la réponse clinique et la clairance plasmatique du médicament. *Gastroentérol Clin Biol*. 1983;7:374-380.
26. Thyss A, Milano G, Renée N, et al. Clinical pharmacokinetic study of 5-FU in continuous 5-day infusions for head and neck cancer. *Cancer Chemother Pharmacol*. 1986;16:64-66.
27. Milano G, Roman P, Khater P, et al. Dose versus pharmacokinetics for predicting tolerance to 5-day continuous infusion of 5-FU. *Int J Cancer*. 1988; 41:537-541.
28. Yoshida T, Araki E, Ligo M, et al. Clinical significance of monitoring serum levels of 5-fluorouracil by continuous infusion in patients with advanced colonic cancer. *Cancer Chemother Pharmacol*. 1990;26:352-354.
29. Trump DL, Egorin MJ, Forrest A, et al. Pharmacokinetic and pharmacodynamic analysis of fluorouracil during 72 hour continuous infusion with and without dipyridamole. *J Clin Oncol*. 1991;9:2027-2035.
30. Fety R, Rolland F, Barberi-Heyob M, et al. Clinical randomized study of 5-FU monitoring versus standard dose in patients with head and neck cancer: preliminary results. *Anticancer Res*. 1994;14:2347-2352.
31. Gamelin EC, Danquechin-Dorval EM, Dumesnil YF, et al. Relationship between 5-fluorouracil dose intensity and therapeutic response in patients with advanced colorectal cancer receiving infusional therapy containing 5-FU. *Cancer*. 1996;77:441-451.
32. Vokes EE, Mick R, Kies MS, et al. Pharmacodynamics of fluorouracil-based induction chemotherapy in advanced head and neck cancer. *J Clin Oncol*. 1996;14:1663-1671.
33. Ychou M, Duffour J, Kramar A, et al. Individual 5-FU dose adaptation in metastatic colorectal cancer: results of a phase II study using a bimonthly pharmacokinetically intensified LV5FU2 regimen. *Cancer Chemother. Pharmacol*. 2003;52(4):282-290.
34. Bertino J, Fleisher M, Beumer JH, et al. Highlights from: 5-fluorouracil drug management pharmacokinetics and pharmacogenomics workshop; Orlando, Florida; January 2007. *Clin Colorectal Cancer*. 2007;6:407-422.
35. Capitain O, Asevoaia A, Boisdrón-Celle M, Soulié P, Morel A, Gamelin E. Influence of pharmacogenetic polymorphisms on 5-fluorouracil and irinotecan efficacy and tolerance in patients treated for advanced colorectal cancer. *Proc Am Soc Clin Oncol Gastrointest Cancers Symp*. [Abstract 429] 2008.
36. Gamelin E, Boisdrón-Celle M, Guerin-Meyer V, Poirier A, Berger V, Morel A. 5-FU dose monitoring and prevention of oxaliplatin-induced neurotoxicity in FOLFOX 4 regimen. Results of a phase II study. *Proc Am Soc Clin Oncol Gastrointest Cancers Symp*. [Abstract 431] 2008.
37. Capitain O, Boisdrón-Celle M, Poirier A-L, Abadie-Lacourtoisie S, Morel A, Gamelin E. The influence of fluorouracil outcome parameters on tolerance and efficacy in patients with advanced colorectal cancer. *Pharmacogenomics J*. 2008;8:256-267. Lee JJ, Beumer JH, Chu E. Therapeutic drug monitoring of 5-fluorouracil. *Cancer chemother and pharmacol*. 2016;78:447-64.
38. Escoriza J, Aldaz A, Calvo E, Giraldez J. Simple and sensitive determination of 5-fluorouracil in plasma by high-performance liquid chromatography. Application to clinical pharmacokinetic studies. *J Chromatogr B Biomed Sci Appl*. 1999;736:97-102.
39. Beumer JH, Chu E, Allegra C, et al. Therapeutic Drug Monitoring in Oncology: International Association of Therapeutic Drug Monitoring and Clinical Toxicology Recommendations for 5-Fluorouracil Therapy. *Clin Pharmacol Ther*. 2019;105(3):598-613.
40. Di Paolo A, Danesi R, Ciofi L, et al. Improved analysis of 5-Fluorouracil and 5,6-dihydro-5-Fluorouracil by HPLC with diode array detection for determination of cellular dihydropyrimidine dehydrogenase activity and pharmacokinetic profiling. *Ther Drug Monit*. 2005;27:362-8.
41. Remaud G, Boisdrón-Celle M, Hameline C, Morel A, Gamelin E. An accurate dihydrouracil/uracil determination using improved high performance liquid chromatography method for preventing fluoropyrimidines-related toxicity in clinical practice. *J Chromatogr B Analyt Technol Biomed Life Sci*. 2005;823:98-107.
42. Kosovec JE, Egorin MJ, Gjurich S, Beumer JH. Quantitation of 5-fluorouracil (5-FU) in human plasma by liquid chromatography/electrospray ionization tandem mass spectrometry. *Rapid Commun Mass Spectrom*. 2008;22:224-230.
43. McNally AJ, Goc-Szkutnicka K, Li Z, Pilcher I, Polakowski S, Salamone SJ. An OnLine Immunoassay for LSD: Comparison with GC-MS and the Abuscreen RIA. *J Anal Toxicol*. 1996;20:404-408.
44. Li Z, Goc-Szkutnicka K, McNally AJ, et al. Synthesis of New d-Propoxyphene Derivatives and the Development of a Microparticle-Based Immunoassay for the Detection of Propoxyphene and Norpropoxyphene. *Bioconj Chem*. 1997; 8:896-905.