

# MyCare Psychiatry Total Aripiprazole Assay Kit

## INDICATIONS FOR USE

The MyCare Psychiatry Total Aripiprazole Assay Kit is intended for the *in vitro* quantitative measurement of total aripiprazole (aripiprazole plus dehydroaripiprazole) in human serum using automated clinical chemistry analyzers. Measurements obtained are used for monitoring patient adherence to aripiprazole therapy to help ensure appropriate treatment.

## SUMMARY AND EXPLANATION OF THE TEST

Aripiprazole (7-[4-[4-(2,3-dichlorophenyl)-1-piperazinyl]butoxy]-3,4-di-hydrocarbostyryl) is a quinolone derivative atypical anti-psychotic agent. It has partial agonistic activity at dopamine D2 receptors and serotonin 5-HT1A receptors and potent antagonistic activity on serotonin 5-HT2A receptors.<sup>1,2</sup> The oral medication is indicated for the treatment of schizophrenia, acute treatment of manic and mixed episodes associated with bipolar disorder, adjunctive treatment of major depressive disorder, irritability associated with autistic disorder, and Tourette's disorder. The injectable is indicated for agitation associated with schizophrenia or bipolar mania. The major metabolite of aripiprazole, dehydroaripiprazole, is also pharmaceutically active.<sup>1</sup> The therapeutic effect of aripiprazole is due to the total exposure to both aripiprazole and the active metabolite (dehydroaripiprazole).<sup>3</sup> The total aripiprazole assay measures the total active aripiprazole in patient serum: aripiprazole plus dehydroaripiprazole.

Nonadherence to medication is well known for patients with severe mental illness.<sup>4</sup> While adherence to medication is critical to successful treatment outcomes, adherence is also least likely to be accurately assessed.<sup>5,6</sup> Measurement of total aripiprazole provides clinicians with objective evidence of concentrations that may be related to patient adherence.<sup>7</sup>

The total aripiprazole assay is a homogenous two reagent nanoparticle agglutination assay used for detection of total aripiprazole in human serum. It is based on competition between drug and drug-conjugates for binding to drug specific antibodies covalently bound to nanoparticles. The extent of particle aggregation can be followed spectrophotometrically on clinical chemistry analyzers.

## REAGENTS

|  |  |                   |
|--|--|-------------------|
| MyCare Psychiatry Total Aripiprazole Assay Kit | REF ARI-RGT  | Quantity x Volume |
| Reagent 1                                      | R1<br>Reaction buffer that contains drug-conjugate, protein and buffer                                     | 1 x 10.0 mL       |
| Reagent 2                                      | R2<br>Nanoparticle reagent that contains monoclonal antibody bound to nanoparticles in a buffered solution | 1 x 5.0 mL        |

## WARNINGS AND PRECAUTIONS

- For In Vitro Diagnostic Use Only.
- For diagnostic purposes, the results should always be assessed with the patient's medical history, clinical examination and other findings.
- Exercise normal precautions required for handling all laboratory reagents.
- Follow reagent handling instructions. Improper mixing of reagents can affect assay performance.
- All components of the total aripiprazole assay contain less than 0.1% sodium azide. Avoid contact with skin and mucous membranes. Flush affected areas with copious amounts of water. Seek immediate medical attention if reagents are ingested or come into contact with eyes. When disposing of such reagents, always flush with large amounts of water to prevent accumulation of azide.

## REAGENT HANDLING

The total aripiprazole assay reagents are ready to use. Mix the reagents (R1 and R2) by gently inverting five times, avoiding the formation of bubbles then place them on the analyzer.

Mix the reagents (R1 and R2) before pouring them into any analyzer-specific (secondary) reagent carrier. Before placing analyzer-specific (secondary) reagent carriers on the analyzer, mix the reagents (R1 and R2) by gently inverting five times, avoiding the formation of bubbles.

## STORAGE AND STABILITY

Store reagents refrigerated at 2 - 8°C. Do not freeze.

When stored and handled as directed, unopened reagents are stable until the expiration date on the label. Improper storage of reagents can affect assay performance.

## SPECIMEN COLLECTION AND HANDLING

Serum is required. Trough or  $C_{min}$  samples at steady state have been recommended for testing antipsychotics.<sup>6</sup> After two weeks of treatment on the same dose, collect samples before the next dose.<sup>8</sup> For long lasting injectables collect the sample before the next dose.<sup>7</sup>

Prepare serum within 3 days of blood collection. Blood and serum samples may be stored at room temperature or 2 - 8°C. Store serum for up to 7 days before measuring. Freeze ( $\leq -20^{\circ}\text{C}$ ) for longer storage. Ensure the sample is thawed and thoroughly mixed before measuring. Avoid repeated freezing and thawing of samples.

## PROCEDURE

### **Materials Provided:**

**REF** ARI-RGT – MyCare Psychiatry Total Aripiprazole Assay Kit

### **Materials Required – Provided Separately:**

**REF** MCP2-CAL – MyCare Psychiatry Calibrator Kit 2

**REF** MCP2-CON - MyCare Psychiatry Control Kit 2

### **Instruments**

Reagents may need to be transferred to analyzer-specific reagent containers.

The performance of applications not validated by Saladax Biomedical, Inc. is not warranted and must be user defined.

### **Assay**

To run the assay, see the instrument specific application sheet and appropriate analyzer operator's manual.

### **Calibration**

Perform a full calibration using the six calibrators in the Calibrator Kit 2. Verify the calibration by testing the low, medium, and high controls in the Control Kit 2.

**Calibration Frequency** - Calibration is recommended:

- After a reagent (kit) lot change,
- After performance of major instrument maintenance,
- As required following quality control procedures.

### **Quality Control (QC)**

Each laboratory should establish its own QC procedures for the total aripiprazole assay kit. All quality control requirements and testing should be performed in accordance with local, state and/or federal regulations or accreditation requirements. Good laboratory practice suggests that at least two QC concentrations be tested each day patient samples are measured, and each time calibration is performed. Ensure that the quality control results meet the acceptance criteria before reporting patient results.

## RESULTS

The concentration result is automatically calculated from the non-linear calibration curve by the analyzer. Report results in ng/mL or nmol/L. The conversion factor from ng/mL is  $2.23 \times \text{ng/mL} = 1 \text{ nmol/L}$ .

## LIMITATIONS OF THE PROCEDURE

The total aripiprazole assay has been validated for serum. Do not use serum separator tubes.

As with any assay utilizing mouse antibodies, the possibility exists for interference by human anti-mouse antibodies (HAMA) in the sample. Samples containing such antibodies can potentially produce erroneous total aripiprazole results, which are inconsistent with the patient's clinical profile.

For samples containing 150 and 500 ng/mL total aripiprazole, 50 ng/mL of cariprazine caused assay biases of 164% and 71% respectively. Elevated levels of aripiprazole may be seen in patients administered cariprazine. For samples containing 150 and 500 ng/mL total aripiprazole, 42,000 ng/mL of lamotrigine (3X the therapeutic level) caused assay biases of 40%. Elevated levels of aripiprazole may be seen in patients administered lamotrigine.

## EXPECTED VALUES

The therapeutic range for total aripiprazole in serum is not fully established. A therapeutic range from 150 to 500 ng/mL has been proposed for aripiprazole plus dehydroaripiprazole.<sup>7</sup> Measured concentrations for adherent patients at steady state are expected to be in the measuring range of the assay. Therapeutic drug monitoring of total aripiprazole has been recommended because of high interpatient variability, unpredictable response, and the importance of adherence for successful therapy.<sup>7</sup> The complexity of the clinical state, individual differences in sensitivity, and co-administered medications may contribute to different requirements for optimal total aripiprazole blood levels. Users should investigate the transferability of the expected values to their own patient population and if necessary determine their own reference range. For diagnostic purposes the test findings should always be assessed in conjunction with the patient's medical history, clinical examinations, and other findings. Clinicians should carefully monitor patients during therapy initiation and dose adjustments. It may be necessary to obtain multiple samples to determine expected variation of optimal (steady state) concentrations for individual patients.

## SPECIFIC PERFORMANCE DATA

Typical performance data for the total aripiprazole assay obtained on a Beckman Coulter AU480 are shown below. Results obtained in individual laboratories may differ from these data.

### **Precision**

Within-laboratory precision and repeatability were verified throughout the measuring range according to CLSI Guideline EP05-A3.<sup>9</sup> Three Control Kit 2 controls (Control 1, 2, 3), two serum pools spiked with both aripiprazole and dehydroaripiprazole to mimic the metabolite ratio found in clinical samples (Serum 1, 2), and two pools of clinical samples (Clinical 1, 2) were tested.

| Sample     | N  | Mean (ng/mL) | Repeatability | Within-Laboratory |
|------------|----|--------------|---------------|-------------------|
|            |    |              | CV            | CV                |
| Control 1  | 80 | 49           | 6.5%          | 8.3%              |
| Control 2  | 80 | 198          | 2.3%          | 4.0%              |
| Control 3  | 80 | 682          | 2.2%          | 3.9%              |
| Serum 1    | 80 | 45           | 6.5%          | 9.5%              |
| Serum 2    | 80 | 959          | 2.6%          | 4.3%              |
| Clinical 1 | 80 | 150          | 3.5%          | 4.1%              |
| Clinical 2 | 80 | 503          | 2.6%          | 4.1%              |

### **Limit of Quantitation (LoQ) and Limit of Detection (LoD)**

The lower limits of quantitation and detection were established using CLSI guideline EP17-A2.<sup>10</sup>

#### **LoQ**

The LoQ was determined with an accuracy goal at the LoQ of  $\leq 35\%$  total error (Westgard model). The LoQ of the total aripiprazole assay is 45 ng/mL.

#### **LoD**

The LoD is the lowest amount of analyte that can be reliably detected ( $\geq 95\%$  of results greater than the limit of blank.). The LoD of the total aripiprazole assay is 22 ng/mL.

### **Result Reporting**

Each laboratory should determine reporting criteria for total aripiprazole concentrations. The following suggestion from CLSI EP17-A2 may be appropriate:<sup>10</sup>

Result  $\leq$  LoB - report "not detected; concentration  $<$  LoD"

LoB  $<$  Result  $<$  LoQ - report "analyte detected; concentration  $<$  LoQ"

Result  $\geq$  LoQ - report the result as measured

## Measurement Range

The measurement range of the total aripiprazole assay is 45 – 1,000 ng/mL.

## Specificity

### Metabolism

Aripiprazole is metabolized in the liver by CYP3A4 and CYP2D6. The major metabolite dehydroaripiprazole also has pharmacological activity.<sup>1,3</sup> At steady state, its concentration is ~40% of the parent drug.<sup>1</sup> The other major metabolite, the acid product of N-dealkylation (OPC-3373) is also present in serum. Another minor metabolite (DCPP) is found at < 20% of the parent drug.

Specificity for the following metabolites was tested in the absence and presence of total aripiprazole at 150, 500, and 1,000 ng/mL.

| Compound   | Tested at (ng/mL) | % Bias |
|--|-------------------|--------|
| 3,4-dihydro-7-(3'carboxy) propoxy-2(1H) quinolinone (OPC-3373) | 475               | 3%     |
| 1-(2,3-dichlorophenyl) piperazine (DCPP)                       | 50                | 6%     |

### Interfering Substances

Testing of interferents was conducted according to CLSI guidelines for interference.<sup>11-13</sup> No significant assay bias was observed from samples with the following endogenous interferents at the given levels:

| Interferent            | Level       |            |
|------------------------|-------------|------------|
| Rheumatoid Factor      | 508 IU/mL   |            |
| Human Serum Albumin    | 10.8 g/dL   | 108 g/L    |
| Human Immunoglobulin G | 12.1 g/dL   | 121 g/L    |
| Icteric Interference   | 43.5 mg/dL  | 744 µmol/L |
| Lipemic Interference   | 614 mg/dL   | 6.9 mmol/L |
| Hemolysate             | 1,050 mg/dL |            |

### Cross-reactivity

Specificity for the following cross-reactants was tested in the absence and presence of total aripiprazole at 150, 500, and 1,000 ng/mL.

Cross-reactivity was tested according to CLSI guidelines for interferences.<sup>11-13</sup> The following compounds did not interfere with the total aripiprazole assay: assay bias was ≤ 13%.

| Compound             | Tested at (ng/mL) | Compound                 | Tested at (ng/mL) |
|----------------------|-------------------|--------------------------|-------------------|
| Acetaminophen        | 200,000           | Acetazolamide            | 60,000            |
| Acetylsalicylic acid | 500,000           | Albuterol                | 1,000             |
| Alendronate sodium   | 1,000             | Alpha - tocopherol       | 130,000           |
| Alprazolam           | 2,000             | Amantadine Hydrochloride | 10,000            |
| Amikacin sulfate     | 144,000           | Amiloride HCl dihydrate  | 500               |
| Amisulpride          | 1,200             | Amitriptyline            | 1,000             |
| Amlodipine besylate  | 100               | S (+)-amphetamine        | 1,000             |
| Amoxapine            | 2,900             | Amoxicillin              | 80,000            |
| L-ascorbic acid      | 60,000            | Asenapine                | 500               |
| Atomoxetine          | 7,900             | Atorvastatin calcium     | 800               |
| Baclofen             | 3,000             | Benztropine              | 600               |
| Betamethasone        | 400               | Biotin                   | 3,600             |
| Biperiden            | 300               | Blonanserin              | 100               |
| Brexipiprazole       | 1,000             | Bromperidol              | 100               |
| Budesonide           | 50                | Bupropion                | 3,000             |

| Compound              | Tested at (ng/mL) | Compound                   | Tested at (ng/mL) |
|-----------------------|-------------------|----------------------------|-------------------|
| Buspirone             | 200               | Caffeine                   | 108,000           |
| Calcium carbonate     | 315,000           | Cannabidiol                | 100               |
| Cannabinol            | 100               | Carbamazepine              | 45,000            |
| L-Carnosine           | 100,000           | Cefalexin                  | 200,000           |
| Celecoxib             | 8,800             | Cetirizine dihydrochloride | 4,400             |
| 8-chloro-theophylline | 3,000             | Chlorpromazine HCl         | 3,300             |
| Cimetidine            | 30,000            | Ciprofloxacin              | 12,000            |
| Citalopram HBr        | 5,500             | Clindamycin                | 51,000            |
| Clonazepam            | 300               | Clotiapine                 | 500               |
| Clotrimazole          | 50                | Clozapine                  | 1,800             |
| Codeine               | 2,000             | Cortisol                   | 300               |
| (-)-cotinine          | 2,000             | Cyclosporin A              | 9,000             |
| Desloratadine         | 600               | Desvenlafaxine             | 800               |
| Dextromethorphan      | 1,000             | Diazepam                   | 30,000            |
| Diphenhydramine HCl   | 6,000             | Divalproex Sodium          | 400,000           |

| Compound                         | Tested at (ng/mL) | Compound                    | Tested at (ng/mL) |
|----------------------------------|-------------------|-----------------------------|-------------------|
| Docosahexaenoic acid ethyl ester | 150,000           | Donepezil                   | 50,000            |
| Doxycycline HCl                  | 35,000            | Droperidol                  | 200               |
| D-Serine                         | 100,000           | Duloxetine                  | 200               |
| Erythromycin                     | 138,000           | Escitalopram                | 200               |
| Estradiol                        | 10                | Eszopiclone                 | 200               |
| Ethanol                          | 10,000,000        | Famotidine                  | 2,500             |
| Fenofibrate                      | 50,000            | Fentanyl                    | 600               |
| Fluoxetine HCl                   | 4,000             | Fluticasone propionate      | 50                |
| Fluvoxamine                      | 2,000             | Folic acid                  | 15                |
| Furosemide                       | 60,000            | Galantamine                 | 200               |
| Gentamicin sulfate               | 30,000            | Glyburide                   | 2,000             |
| Haloperidol                      | 1,000             | Heparin sodium salt         | 50 U/mL           |
| Hydrochlorothiazide              | 6,000             | Hyoscine (Scopolamine HBr)  | 100               |
| Hyperforin (St. John's Wort)     | 200               | Hypericin (St. John's Wort) | 100               |
| Ibuprofen                        | 500,000           | lloperidone                 | 100               |
| Imipramine                       | 700               | Indinavir sulfate           | 400               |
| Lactulose                        | 10,000            | Lamivudine                  | 10,500            |
| Lamotrigine                      | 15,000            | Lansoprazole                | 9,400             |
| Levonorgestrel                   | 100               | Lisinopril dihydrate        | 350               |
| Lithium carbonate                | 250,000           | Lorazepam                   | 1,000             |
| Lovastatin                       | 500               | Loxapine                    | 300               |
| Lurasidone                       | 400               | Meclizine dihydrochloride   | 500               |
| Metformin                        | 40,000            | Methotrimeprazine           | 600               |
| Methylphenidate HCl              | 350               | Metoclopramide HCl          | 500               |
| Metoprolol tartrate              | 5,000             | Metronidazole               | 123,000           |
| Midazolam                        | 3,800             | Milnacipran                 | 10,000            |
| Mirtazapine                      | 900               | Mometasone furoate          | 50                |
| Morphine                         | 7,800             | Naltrexone                  | 200               |
| Naproxen sodium                  | 500,000           | Nateglinide                 | 30,000            |
| Nefazodone HCl                   | 6,000             | Nicotine                    | 1,000             |
| Nicotinic acid                   | 27,900            | Nordiazepam                 | 5,000             |
| Nortriptyline                    | 1,200             | Olanzapine                  | 300               |
| Omeprazole                       | 8,400             | Oxazepam                    | 5,000             |

| Compound                | Tested at (ng/mL) | Compound                    | Tested at (ng/mL) |
|-------------------------|-------------------|-----------------------------|-------------------|
| Oxcarbazepine           | 105,000           | Oxycodone                   | 500               |
| Paliperidone            | 60                | Pantothenic acid            | 1,800             |
| Paroxetine              | 1,200             | Penicillin V                | 42,000            |
| Perazine                | 1,400             | Perlapine                   | 150               |
| Perphenazine            | 100               | Phenobarbital               | 690,000           |
| Phentermine             | 500               | Phenytoin                   | 60,000            |
| Pimozide                | 100               | Pipamperone dihydrochloride | 1,200             |
| Potassium EDTA          | 1000              | Pravastatin sodium          | 300               |
| Prednisolone            | 3,000             | Pregabalin                  | 22,500            |
| Procyclidine            | 1,900             | Promethazine                | 1,200             |
| R,R (-)-pseudoephedrine | 10,000            | S,S (+)-pseudoephedrine     | 10,000            |
| Pyridoxine HCl          | 100               | Quetiapine                  | 2,800             |
| Quinidine               | 15,000            | Raloxifene                  | 50                |
| Ranitidine              | 10,500            | Retinol                     | 4,000             |
| Riboflavin              | 200               | Rifampicin                  | 65,000            |
| Risperidone             | 200               | Rosuvastatin calcium        | 200               |
| Salicylic acid          | 500,000           | Sarcosine                   | 1,500             |
| Sertindole              | 300               | Sertraline hydrochloride    | 1,000             |
| Simvastatin             | 1,700             | Sodium benzoate             | 400,000           |
| Sodium fluoride         | 900               | Spiro lactone               | 600               |
| Sulfamethoxazole        | 400,000           | Sulpiride                   | 50,000            |
| Temazepam               | 5,000             | Terbinafine                 | 9,000             |
| Theophylline            | 60,000            | Thiamine HCl                | 500               |
| Topiramate              | 75,000            | Trazodone HCl               | 14,000            |
| Triamcinolone acetonide | 300               | Triamterene                 | 9,000             |
| Triazolam               | 40                | Valproic acid               | 500,000           |
| Vancomycin HCl          | 120,000           | Varenicline                 | 50                |
| Venlafaxine HCl         | 700               | Vitamin B12                 | 50                |
| Vitamin D2              | 200               | Vitamin K1                  | 50                |
| Warfarin                | 75,000            | Ziprasidone                 | 600               |
| Zolpidem hemitartrate   | 5,000             | Zonisamide                  | 120,000           |
| Zopiclone               | 200               | Zuclopenthixol              | 300               |

### Recovery

The recovery of total aripiprazole was assessed in the 3 controls, two spiked serum pools and two clinical pools measured for the EP05-A3 precision performance study. The percent recovery was determined by dividing the mean measured concentration of each sample by the expected concentration total aripiprazole. All mean recoveries were within 88% to 114%.


### Linearity

The linearity of the total aripiprazole assay was verified according to CLSI guideline EP6-A.<sup>14</sup> Eleven linearity samples covering the measuring range were prepared in human serum spiked with aripiprazole. Deviation from linearity (n=5) was ≤ 10%. The assay was linear across the measuring range from 45 – 1,000 ng/mL.

### Method Comparison

Results of the total aripiprazole assay were compared to a validated LC-MS/MS according to CLSI guideline EP09-A3.<sup>15</sup> Deming regression analysis was performed with 110 patient samples. Results are shown for one lot.














| Deming Regression Statistics<br>Total Aripiprazole Assay vs. LC-MS/MS |          |
|---|----------|
| Slope   | 1.01     |
| Intercept   | 2.56     |
| Correlation Coefficient (R)   | 0.98     |
| N   | 110      |
| Concentration Range (LC-MS/MS)  | 48 - 839 |



## References

- Otsuka America Pharmaceutical I. Abilify (Aripiprazole) Prescribing Information. Product Insert. 2017.
- PubChem Aripiprazole <https://pubchem.ncbi.nlm.nih.gov/compound/60795> accessed March 30, 2017
- Lin SK, Chen CK, Liu YL. Aripiprazole and dehydroaripiprazole plasma concentrations and clinical responses in patients with schizophrenia. *J Clin Psychopharmacol.* 2011;31(6):758-762.
- Velligan DI, Weiden PJ, Sajatovic M, et al. Assessment of adherence problems in patients with serious and persistent mental illness: recommendations from the Expert Consensus Guidelines. *J Psychiatr Pract.* 2010;16(1):34-45.
- Higashi K, Medic G, Littlewood KJ, Diez T, Granstrom O, De Hert M. Medication adherence in schizophrenia: factors influencing adherence and consequences of nonadherence, a systematic literature review. *Ther Adv Psychopharmacol.* 2013;3(4):200-218.
- Haddad PM, Brain C, Scott J. Nonadherence with antipsychotic medication in schizophrenia: challenges and management strategies. *Patient Relat Outcome Meas.* 2014;5:43-62.
- Hiemke C, Bergemann N, Clement HW, et al. Consensus Guidelines for Therapeutic Drug Monitoring in Neuropsychopharmacology: Update 2017. *Pharmacopsychiatry.* 2018;51:9-62.
- Grundmann M, Kacirova I, Urinovska R. Therapeutic drug monitoring of atypical antipsychotic drugs. *Acta Pharm.* 2014;64(4):387-401.
- CLSI. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. CLSI document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute, 2014.
- CLSI. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- CLSI. Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition CLSI document EP7-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
- CLSI. Interference Testing in Clinical Chemistry. 3rd ed. CLSI guideline EP07. Wayne PA: Clinical and Laboratory Standards Institute; 2018.
- CLSI. Supplemental Tables for Interference Testing in Clinical Chemistry. 1st ed. CLSI supplement EP37. Wayne PA: Clinical and Laboratory Standards Institute; 2018.
- NCCLS. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. NCCLS document EP6-A. Wayne, PA: NCCLS; 2003.
- CLSI. Measurement Procedure and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition. CLSI document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.

**SYMBOLS USED**

|  |  |   |   |
|--|--|---|---|
|   | <i>in vitro</i> Diagnostic Device        |        | Consult Instructions for Use                                      |
|   | Catalog Number                           |        | Use By  |
|   | Batch Code                               |        | Temperature Limitation  |
|   | Manufacturer                             | <b>Rx only</b>  | For Prescription Use Only   |
| <br> | Reagent 1<br>Reagent 2                   |  (N) x | Gently invert reagents (R1 and R2) N number of times prior to use |
|   | CE mark                                  |        | Authorized Representative in the European Community               |
|   | Authorized Representative in Switzerland |   |   |



**EMERGO EUROPE**  
Westervoortsedijk 60,  
6827 AT Arnhem  
The Netherlands



Casus Switzerland GmbH  
Hinterbergstrasse 49  
6312 Steinhausen  
Switzerland

Saladax Biomedical, Inc.  
116 Research Dr.  
Bethlehem, PA 18015 USA  
MyCareTests.com

Customer Service  
Telephone: +1 (610) 419-6731  
Fax: +1 (484) 547-0590  
Email: Techsupport@saladax.com

© 2026, Saladax Biomedical, Inc.

MyCare™ is a trademark of Saladax Biomedical, Inc. All other product names and trademarks are property of their respective owners.