



MyCare Psychiatry Assays Clinical Chemistry Analyzers

CANADA PRODUCT CATALOG

| Reagents | Order Code | Description | Shelf Life | Contents |
|-----------------------------------|------------|------------------|------------|---|
| Total Aripiprazole | ARI-RGT | Reagent Kit | 10 Months | 1 vial of Reagent 1 (10mL) and Reagent 2 (5mL); ~100 Tests* |
| Clozapine | CLZ-RGT | Reagent Kit | 10 Months | 1 vial of Reagent 1 (10mL) and Reagent 2 (5mL); ~100 Tests* |
| Olanzapine | OLZ-RGT | Reagent Kit | 10 Months | 1 vial of Reagent 1 (10mL) and Reagent 2 (5mL); ~100 Tests* |
| Paliperidone | PAL-RGT | Reagent Kit | 10 Months | 1 vial of Reagent 1 (10mL) and Reagent 2 (5mL); ~100 Tests* |
| Total Risperidone | RSP-RGT | Reagent Kit | 10 Months | 1 vial of Reagent 1 (10mL) and Reagent 2 (5mL); ~100 Tests* |
| Quetiapine | QTP-RGT | Reagent Kit | 10 Months | 1 vial of Reagent 1 (10mL) and Reagent 2 (5mL); ~100 Tests* |
| Calibrators & Controls | MCP1-CAL | Calibrator Kit 1 | 10 Months | 1 set of 6 calibrators (A – F); 3mL per vial; For calibration of total risperidone and paliperidone assays |
| | MCP1-CON | Control Kit 1 | 10 Months | 2 sets of 3 control levels (Low, Medium, High); 3mL per vial; For use with total risperidone and paliperidone assays |
| | MCP2-CAL | Calibrator Kit 2 | 10 Months | 1 set of 6 calibrators (A – F); 3mL per vial; For use with clozapine, total aripiprazole, quetiapine, and olanzapine assays |
| | MCP2-CON | Control Kit 2 | 10 Months | 2 sets of 3 control levels (Low, Medium, High); 3mL per vial; For use with clozapine, total aripiprazole, quetiapine, and olanzapine assays |

*Analyzer Dependent